

Incubator Hood TH 15



Description

The Incubator Hood TH 15 has a temperature range from +5°C above ambient to +50°C and can be combined with all compact shakers KS 15 and Microplate Shakers TiMix 5 or used as independent unit. The height is sufficient for 1000 ml Erlenmeyer flasks.

Basic Equipment

Coated metal housing with 3 windows, the front door opens upwards

Safety cut-off of the shaker if combined with a Bühler Compact Shaker

Standard temperature controller 5°C above ambient to +50°C

Option

[Extended temperature range up to +60°C](#)

Details

Temperature range:	5°C above ambient to +50°C
Temperature accuracy:	+ / - 1% of the set value
Inner dimensions (w x d x h):	90 x 410 x 310 mm (height sufficient for 1000 ml Erlenmeyer flasks)
Order no.:	6161 000
Accessory:	Einhängetablar TH 15

Technical specifications

Temperature sensor:	PT 100
Air circulation:	240 m ³ /h
Electrical supply:	230 V , 50 / 60 Hz (115 V version on request)
Dimensions (w x d x h):	510 x 450 x 455 mm
Weight:	21 kg