

Infrared Heater IRB 1



Description

The compact laboratory heater IRB 1 with a heating area of 60 x 60 mm is optimally suited for fast heating tasks without open flame up to max. 700°C.

Basic equipment

Basic device with support cross

Details

Max. temperature:	700°C
Heating power:	250 W
Heating area:	60 x 60 mm
Max. ambient temperature:	30°C
Order no.:	6060 000
Accessory:	Leistungsregler LR 1 Stativhalterung

Technical specifications

Electrical supply:	230 V, 50 / 60 Hz (115 V version on request)
Dimensions (w x d x h):	100 x 105 x 105 mm
Weight:	0,7 kg