

Infrared Heater IRB 2



Description

The compact laboratory heater IRB 2 with power control and a heating area of 100 x 100 mm is optimally suited for fast heating tasks without open flame up to max. 900°C.

Basic equipment

Basic device incl. integrated power control, support cross

Details

Max. temperature:	900°C
Heating power:	800 W
Heating area:	100 x 100 mm
Power control:	stepless control between 0- 100%
Max. ambient temperature:	30°C
Order number:	6063 000
Accessory:	Stativhalterung

Technical specifications

Electrical supply:	230 V, 50 / 60 Hz (115 V version on request)
Dimensions (w x d x h):	150 x 150 x 170 mm
Weight:	2,5 kg