

## Dual-Action Shaker KL 2



### Description

Depending on the task the shaking motion of the Dual Action Shaker KL 2 can be switched from orbital to reciprocating.

#### Basic Equipment

Basic device incl. rack system with 3 clamping strips and rubber mat

#### Loading Capacity

Erlenmeyer flasks 250 ml 9 pieces

Separating funnels 250 ml 1 pieces

You will find further loading possibilities [here](#)

## Details

<b>Motion:</b>	can be switched from orbital to reciprocating (horizontal)
<b>Shaking platform:</b>	300 x 300 mm
<b>Max. load:</b>	2 kg
<b>Shaking speed:</b>	30 - 420 rpm, steplessly adjustable
<b>Stroke:</b>	8 mm
<b>Timer:</b>	0 - 120 min / continuous
<b>Order No.:</b>	6115 000
<b>Accessory:</b>	Spannleiste v 0050 477 Spannleiste h 0050 118 Klammerleiste 0050 206 Federleiste 0050 207 Reagenzglasgestell mit 44 Bohrungen à 14 mm Ø Reagenzglasgestell mit 44 Bohrungen à 16 mm Ø Reagenzglasgestell mit 44 Bohrungen à 18 mm Ø Reagenzglasgestell mit 14 Bohrungen à 30 mm Ø

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 Hz please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 4 W
<b>Ambient temperature:</b>	5°C bis 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	355 x 455 x 195 mm
<b>Weight:</b>	13,5 kg