


# Microplate Shaker TiMix 2



## Description

The Microplate Shaker TiMix 2 especially for small samples in microplates, deepwell-plates or other plates in microplate format.

### Basic Equipment

Basic device incl. basic platform with holders for 4 microplates

### Loading Capacity

4 microplates

(deepwell-plates, or other plates in microplate format)

You will find further loading possibilities [here](#)

## Details

<b>Motion:</b>	exact orbital motion
<b>Shaking platform:</b>	200 x 295 mm for 4 microplates
<b>Max. load:</b>	2 kg
<b>Shaking speed:</b>	100 - 1400 min <sup>-1</sup> , steplessly adjustable
<b>Stroke:</b>	3 mm
<b>Timer:</b>	0 - 120 min / continuous
<b>Order No.:</b>	6110 000
<b>Accessory:</b>	Schwingaufsatz KM Spannleiste KM

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 Hz please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 6 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	220 x 410 x 110 mm
<b>Weight:</b>	8,5kg