


Three-Dimensional Shaker TL 10



Description

The Three-Dimensional Shaker TL 10 works with slow, three-dimensional motion. The angle of tilt is adjustable between 0° and 9° and the shaking speed between 1 and max. 23 rpm. The shaking platform of 340 x 340 mm can be loaded with max. 10 kg.

Basic Equipment

Basic device incl. basic tray and rubber mat

Details

Motion: three-dimensional orbital motion

Shaking platform: 340 x 340 mm
Max. load: 10 kg
Shaking speed: 1 - 23 rpm, steplessly adjustable
Angle of tilt: 0° to 9° einstellbar
Timer: 0 - 120 min / continuous
Order No.: 6117 000
Accessory: Mehretagen-Aufsatz WS / TL

Technical specifications

Electrical supply: 230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection: IP 21
Heat emission: approx 6 W
Ambient temperature:: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 355 x 450 x 205 mm
Weight: 18,5 kg