

Tilt Shaker WS 10



Description

The Tilt Shaker WS 10 works with slow, linear tilting motion and a fix angle of tilt of 5° with a shaking speed of 1 - 23 rpm. The shaking platform of 340 x 340 mm can be loaded with a weight of 10 kg.

Basic Equipment

Basic device incl. basic tray and rubber mat

Details

Motion:	linear tilting motion
Shaking platform:	340 x 340 mm
Max. load:	10 kg
Shaking speed:	1 - 23 rpm, steplessly adjustable
Angle of tilt:	5° fixed
Timer:	0 - 120 min / continuous
Order No.:	6114 000
Accessory:	Mehretagen-Aufsatz WS / TL

Technical specifications

Electrical supply:	115 - 230 V, 50 / 60 Hz
Enclosure protection:	IP 21
Heat emission:	approx. 3 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	355 x 450 x 205 mm
Weight:	18,5 kg