

# Tube Roller RM 2



## Description

Gentle tilting and rotating motion for uniform mixing

### Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

### Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

## Details

<b>Motion:</b>	rotating and tilting motion
<b>Shaking platform:</b>	200 x 300 mm

**Max. load:** 2 kg  
**Shaking speed:** 5 - 80 rpm, in steps of 5  
**Angle of tilt:** 3° fixed  
**Stroke at end of rolls:** 15 mm  
**Dimensions of rolls:** 300 x Ø 32 mm  
**Number of rolls:** 6 pieces  
**Runtime:** 0 - 99 h 59 min / continuous  
**Order No.:** 6138 000  
**Accessory:** Überkopf Aufsatz RM 2

## Technical specifications

**Electrical supply:** 100 - 240 V, 50 / 60 Hz  
**Enclosure protection:** IP 21  
**Heat emission:** max. 0,5 - 2 W  
**Ambient temperature:** 5°C to 45°C  
**Relative humidity:** ~ 85 % (non-condensing)  
**Dimensions (w x d x h):** 244 x 465 x 115 mm  
**Weight:** 7,3 kg