

# Beer Bottle Shaker SM 30 B control



## Description

The laboratory shaker SM 30 B control with a horizontal shaking motion up to 240 rpm and a stroke of 50 mm is ideally suited for Gushing Tests. It also has a touch display with programmable control. Up to 6 horizontal bottles (0,33 l or 0,5 l) can be fixed on the shaker with the suitable rack system. With a loading capacity of max. 30 kg it can be supplied with a two- or three-storey top frame.

### Basic Equipment

Basic device with shaking plate and rubber mat, **without rack system**  
With touch display and programmable control

### Loading Capacity

Up to 6 horizontal bottles (0,33 or 0,5l) per storey

1-storey, 2-storey or 3-storey rack system can be selected as accessories

[Loading details](#) against shaking speed.

## Options

[USB 2.0 interface](#)

[RS 232 interface](#)

[Control interface](#)

[Safety Package](#)

[Defined stop](#)

## Details

<b>Motion:</b>	reciprocating (horizontal)
<b>Shaking platform:</b>	560 x 400 mm
<b>Max. load:</b>	30 kg
<b>Shaking speed:</b>	up to 240 rpm, in steps of 5
<b>Stroke:</b>	50 mm for Gushing Test (No. 0052 216)
<b>Runtime:</b>	programmable / continuous
<b>Order No.:</b>	6102 000
<b>Accessory:</b>	1-etagiger Aufsatz 2-etagiger Aufsatz 3-etagiger Aufsatz

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 Hz please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 8 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	680 x 610 x 160 mm
<b>Weight:</b>	33 kg