

Barrel Shaker VKS 75 F control



Description

The big heavy-duty shaker for barrels VKS 75 F control with touch display and programmable control has a horizontal shaking motion. With a shaking speed up to 80 rpm, a selected stroke between 120 - 150 mm, a loading capacity up to 75 kg and a respective rack system with 3 fixing elements the shaker is optimally suited for horizontal or vertical barrels up to 80 l capacity. **(floor installation only!)**

Basic Equipment

Basic device **without shaking plate, without rack system**

With touch display and programmable control

Loading Capacity

Horizontal or vertical barrels up to 80 l capacity

[Loading details](#) against shaking speed.

Options:

[USB 2.0 interface](#)

[RS 232 interface](#)

[Control interface](#)

[Safety Package](#)

Details

Motion:	reciprocating (horizontal)
Shaking platform:	760 x 600 mm
Max. load:	75 kg
Shaking speed:	up to 80 rpm, in steps of 5
Stroke:	120 -150 mm in 10 mm steps
Timer:	programmable / continuous
Order no.:	6212 000
Accessory:	Schwingplatte VKS (1-etagig) Sonderaufsatz für Fässer Sonderaufsatz für liegende Fässer

Technical specifications

Electrical supply:	100 - 240 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Heat emission:	max. 15 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	1050 x 835 x 250 mm
Weight:	110 kg
Total weight (Shaker + shaking plate + Rack system):	132 kg