


Champagne Bottle Shaker



Description

The big heavy-duty shaker VKS 75 B control with touch display and programmable control has a horizontal shaking motion. With a shaking speed up to 80 rpm, a stroke of 130 mm, a loading capacity up to 75 kg and an adapted rack system the shaker is optimally suited for up to 36 horizontal champagne bottles or bottles of similar size. **(floor installation only!)**

Basic Equipment

Basic device without shaking plate, without rack system.
With touch display and programmable control

Loading Capacity

Up to 36 horizontal champagne bottles.

Options:

[USB 2.0 interface](#)

[RS 232 interface](#)

[Control interface](#)

[Safety Package](#)

Details

Motion:	reciprocating (horizontal)
Shaking platform:	760 x 600 mm
Max. load:	75 kg
Shaking Speed:	up to 80 rpm, in steps of 5
Stroke:	130 mm
Timer:	programmable / continuous
Order no.:	6112 S01
Accessory:	Aufsatz für Sektflaschen

Technical specifications

Electrical supply:	100 - 240 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Heat emission:	max. 15 W
Umgebungstemperatur:	5°C bis 50°C
Ambient temperature:	~ 85 %
Dimensions shaker (w x d x h)::	1050 x 835 x 250 mm
Weight:	110 kg
Total weight (Shaker + Rack system):	135 kg