

## Loading Capacity of Universal trays

with spring clamps (stainless steel) - maximum numbers

Size	KM Mini 200 x 295 mm	KM 356 x 300 mm	KS 430 x 300 mm	SM 560 x 400 mm	VKS 760 x 600 mm
10 ml	23	60	70	118	279
25 ml	15	30	35	64	132
50 ml	15	30	35	63	132
100 ml	7	16	20	31	63
250 ml	4	12	15	24	46
500 ml	2	9	12	14	34
1000 ml	1	5	5	11	20
2000 ml	1	2	2	6	12
3000 ml	-	-	2	4	8
5000 ml	-	-	1	2	6

with test tube racks (with hinged foot) - maximum numbers

	KM Mini 200 x 295 mm	KM 356 x 300 mm	KS 430 x 300 mm	SM 560 x 400 mm	VKS 760 x 600 mm
	1	3	3	5	8

## Loading Capacity of rack systems with clamping strips

with Erlenmeyer flasks - maximum numbers



Size	KM 200 x 295 mm	KL 300 x 300 mm	KS 400 x 275 mm	SM 560 x 400 mm	VKS 760 x 600 mm
10 ml	20 <sup>1)</sup>				
25 ml	16 <sup>1)</sup>				
50 ml	12 <sup>1)</sup>	20 <sup>5)</sup>	25 <sup>7)</sup>	49 <sup>8)</sup>	104 <sup>10)</sup>
100 ml	9 <sup>2)</sup>	12 <sup>6)</sup>	20 <sup>7)</sup>	42 <sup>8)</sup>	82 <sup>11)</sup>
250 ml	6 <sup>2)</sup>	9 <sup>6)</sup>	12 <sup>5)</sup>	20 <sup>9)</sup>	45
500 ml	2 <sup>3)</sup>	4	6 <sup>6)</sup>	12	32
1000 ml	2 <sup>3)</sup>	2	4	6	20
2000 ml	1 <sup>4)</sup>	1	2	6	12
3000 ml	1 <sup>4)</sup>	1	1	4	9
5000 ml	-	1	1	2	6

<sup>1)</sup> with 5 clamping strips h (0052 115) additional <sup>2)</sup> with 4 clamping strips h (0052 115) additional <sup>3)</sup> with 3 clamping strips h (0052 115) additional <sup>4)</sup> with 2 clamping strips h (0052 115) additional <sup>5)</sup> with 2 clamping strips h (0050 118) additional <sup>6)</sup> with 1 clamping strip h (0050 118) additional <sup>7)</sup> with 3 clamping strips h (0050 118) additional <sup>8)</sup> with 3 clamping strips h (0050 400) additional <sup>9)</sup> with 1 clamping strip h (0050 400) additional <sup>10)</sup> with 5 clamping strips h (0050 387) additional <sup>11)</sup> with 1 clamping strip h (0050 387) additional

with test tube racks - maximum numbers

	KM 200 x 295 mm	KL 300 x 300 mm	KS 400 x 300 mm	SM 560 x 400 mm	VKS 760 x 600 mm
	1	1	2	4	8

with separating funnels - maximum numbers

Size	KM 200 x 295 mm	KL 300 x 300 mm	KS 400 x 275 mm	SM 560 x 400 mm	VKS 760 x 600 mm
100 ml		4 <sup>12)</sup>	4 <sup>12)</sup>	10	14
250 ml		1 <sup>12)</sup>	1 <sup>12)</sup>	6	10
500 ml			1 <sup>12)</sup>	2	6 - 8 <sup>13)</sup>
1000 ml			1 <sup>12)</sup>	2	6
2000 ml			1 <sup>12)</sup>	1 - 2 <sup>13)</sup>	2 - 4 <sup>13)</sup>

<sup>12)</sup> with 1 clamp strip (0050 206) and 1 spring strip (0050 207) additional <sup>13)</sup> Number of pieces depends on shape and dimensions of the separating funnels

Differing loading possibilities on request.